



NATIONAL INSTITUTE OF MENTAL HEALTH & NEURO SCIENCES
INSTITUTE OF NATIONAL IMPORTANCE
BENGALURU - 560 029 INDIA

Date: 11.04.2019

MERIT LIST FOR POSTDOCTORAL FELLOWSHIP IN TRANSFUSION MEDICINE

COURSE

The selection procedure for admission to this course will be followed as indicated in the Prospectus 2019-20.

Individual communication will be sent by email **ONLY** to the provisionally selected candidates on or before **1st May 2019**.

The seats will be allotted in the order of merit only.

In the case of two or more candidates obtaining equal marks in the Online Entrance test, their inter se merit will be determined as indicated in Page 73 of the Prospectus 2019-20 and candidate ranked higher will only be given the seat.

All provisionally selected candidates willing to take provisional admission to various PG courses should upload documents/certificates in the willingness window provided on NIMHANS website between 20th April 2019 to 30th April 2019 failing which, the seat will be allotted to the next candidate in the order of merit.

Candidates selected for multiple courses should send an email giving the preference of one course which he/she wants to join at NIMHANS latest by 25th April 2019 failing which, their candidature will not be considered for all the courses selected.

No individual calls or emails regarding admission to the seats offered will be entertained by NIMHANS.

Candidates who have already completed one Postdoctoral Fellowship course at NIMHANS shall not be eligible/allowed to apply/join for any other Postdoctoral Fellowship course at NIMHANS as indicated in the Prospectus 2019-20.

Please note: The selection is provisional for candidates shortlisted only and admission is subject to verification of results, original documents, tie-breaking and all other criteria indicated in the Prospectus 2019-20 on 4th June 2019.

<i>Prov. Reg. No.</i>	<i>Roll Number</i>	<i>Total Score</i>
19201007196	71323	68.75
19201004001	71294	66.25
19201002191	76658	58.25
19201006243	71310	52
19201002175	71283	38.5

Sd/-
DIRECTOR